# DNBelab\_C\_Series\_HT\_scRNA-analysis-software 2.0流程、安装与使用

# CONTENTS 目录

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
# DNBelab\_C4\_scRNA\_V2技术原理

## Droplet-based indexing

## Deconvolution based on barcode and M280 Index

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









	reads	time	average
	QC pass		
1-1-cDNA-20220523	586,705,374	7:11	
2-1-cDNA-20220523	489,356,785	6:17	
3-1-cDNA-20220523	524,144,512	6:52	7小时
6-1-cDNA-20220523	532,351,518	7:12	
8-1-cDNA-20220523	462,329,125	6:25	
8-2-cDNA-20220523	450,914,187	7:13	

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



Reads mapped to genome (Map Quality $\geq 0$ )	96.29%
Reads mapped to exonic regions	83.8%
Reads mapped to intronic regions	3.2%
Reads mapped antisense to gene	5.9%
Reads mapped to intergenic regions	13.0%
Include introns	True

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Including intronic reads, for both cellular and nuclei samples, could lead to higher sensitivity (higher UMI counts, more genes per cell) and less wasted sequencing.

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
Cellranger方法：高阈值+emptydrops（默认）

EmptyDrops: distinguishing cells from empty droplets in droplet-based single-cell RNA sequencing data

Aaron T. L. Lun<sup>1,†</sup>, Samantha Riesenfeld<sup>2,†</sup>, Tallulah Andrews<sup>3,†</sup>, The Phuong Dao<sup>4,†</sup>, Tomas Gomes<sup>3,†</sup>, participants in the 1<sup>st</sup> Human Cell Atlas Jamboree and John C. Marioni<sup>1,3,5\*</sup>

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R

```
#' Read feature count matrix generated by `PISA count`.
#'
#' This function will read Matrix Market files from a directory which generated by `PISA count`.
#' @import Matrix
#' @param mex_dir Feature count outdir generated by `PISA count`.
#' @return Returns a sparse matrix of feature counts or a list of spliced, unspliced, and
#'         spanning reads sparse matrix.
#'
#' @export
ReadPISA <- function(mex_dir=NULL,
                     barcode.path = NULL,
                     feature.path = NULL,
                     matrix.path=NULL,
                     use_10X=FALSE) {
  if (is.null(mex_dir) && is.null(barcode.path) && is.null(feature.path) &&
      is.null(matrix.path)) {
    stop("No matrix set.")
  }
  if (!is.null(mex_dir) && !file.exists(mex_dir) ) {
    stop(paste0(mex_dir, " does not exist."))
  }
  if (is.null(barcode.path) && is.null(feature.path) && is.null(matrix.path)) {
    barcode.path <- paste0(mex_dir, "/barcodes.tsv.gz")
    feature.path <- paste0(mex_dir, "/features.tsv.gz")
    matrix.path <- paste0(mex_dir, "/matrix.mtx.gz")
  }
  spliced.path <- paste0(mex_dir, "/spliced.mtx.gz")
  unspliced.path <- paste0(mex_dir, "/unspliced.mtx.gz")
  spanning.path <- paste0(mex_dir, "/spanning.mtx.gz")

  if (!file.exists(barcode.path) || !file.exists(feature.path)) {
    stop(paste0("No expression file found at ", mex_dir))
  }
}
```

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2

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GitHub官网:

[https://github.com/MGI-tech-bioinformatics/DNBelab\\_C\\_Series\\_HT\\_scRNA-analysis-software](https://github.com/MGI-tech-bioinformatics/DNBelab_C_Series_HT_scRNA-analysis-software)

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## DNBelab\_C4\_scRNA\_V2流程安装

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运行 `git clone` [https://github.com/MGI-tech-](https://github.com/MGI-tech-bioinformatics/DNBelab_C_Series_HT_scRNA-analysis-software.git)

[bioinformatics/DNBelab\\_C\\_Series\\_HT\\_scRNA-analysis-software.git](https://github.com/MGI-tech-bioinformatics/DNBelab_C_Series_HT_scRNA-analysis-software.git)

`chmod 755 -R DNBelab_C_Series_HT_scRNA-analysis-software`更改权限

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运行 `sh Miniconda3-py37_4.9.2-Linux-x86_64.sh`

```
Do you accept the license terms? [yes|no]
[no] >>> yes
```

```
Preparing transaction: done
Executing transaction: done
installation finished.
Do you wish the installer to initialize Miniconda3
by running conda init? [yes|no]
[no] >>> yes
```

```
Do you accept the license terms? [yes|no]
[no] >>> yes
```

```
Miniconda3 will now be installed into this location:
/home/mgi/miniconda3
```

- Press ENTER to confirm the location
- Press CTRL-C to abort the installation
- Or specify a different location below

```
[/home/mgi/miniconda3] >>> /home/mgi/miniconda3
```

==> For changes to take effect, close and re-open your current shell. <==

If you'd prefer that conda's base environment not be activated on startup, set the `auto_activate_base` parameter to false:

```
conda config --set auto_activate_base false
```

Thank you for installing Miniconda3!

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```
wget https://github.com/broadinstitute/cromwell/releases/download/35/cromwell-35.jar
```

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```
DNBC4tools mkref --action mkgtf --ingtf gene.gtf --outgtf gene.filter.gtf \
--attribute gene_type:protein_coding \
    gene_type:lncRNA \
    gene_type:IG_C_gene \
    gene_type:IG_D_gene \
    gene_type:IG_J_gene \
    gene_type:IG_LV_gene \
    gene_type:IG_V_gene \
    gene_type:IG_V_pseudogene \
    gene_type:IG_J_pseudogene \
    gene_type:IG_C_pseudogene \
    gene_type:TR_C_gene \
    gene_type:TR_D_gene \
    gene_type:TR_J_gene \
    gene_type:TR_V_gene \
    gene_type:TR_V_pseudogene \
    gene_type:TR_J_pseudogene
```

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- 修改run.sh路径, 并填写config.json
- 运行run.sh

run.sh:

```
export PATH=/Local/path/miniconda3/envs/DNBC4tools/bin:$PATH
export LD_LIBRARY_PATH=/Local/path/miniconda3/envs/DNBC4tools/lib:$LD_LIBRARY_PATH
java -jar /Local/path/pipeline/workflows/cromwell-35.jar run -i config.json /Local/path/pipeline/wdl/DNBC4_scRNA.wdl
```

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config.json

```
{
  "main.Outdir": "/Local/path/to/outdir/pbmc_demo",
  "main.SampleName": "pbmc_demo",
  "main.cDNA_Fastq1": "/Local/path/data/cDNA/L01_read_1.fq.gz,/Local/path/data/cDNA/L02_read_1.fq.gz",
  "main.cDNA_Fastq2": "/Local/path/data/cDNA/L01_read_2.fq.gz,/Local/path/data/cDNA/L02_read_2.fq.gz",
  "main.Oligo_Fastq1": "/Local/path/data/oligo/L01_read_1.fq.gz,/Local/path/data/oligo/L02_read_1.fq.gz",
  "main.Oligo_Fastq2": "/Local/path/data/oligo/L01_read_2.fq.gz,/Local/path/data/oligo/L02_read_2.fq.gz",
  "main.BeadsBarcode": "/Local/path/to/pipeline/DNBC4tools/config/DNBelabC4_scrRNA_beads_readStructure.json",
  "main.OligoBarcode": "/Local/path/to/pipeline/DNBC4tools/config/DNBelabC4_scrRNA_oligo_readStructure.json",
  "main.Root": "/Local/path/to/pipeline",
  "main.Refdir": "/Local/path/to/GRCh38/star_index",
  "main.Gtf": "/Local/path/to/GRCh38/genes.filter.gtf",
  "main.Oligo_type8": "/Local/path/to/pipeline/DNBC4tools/config/oligo_type8.txt",
  "main.Species": "Human",
  "main.expectCellNum": 3000,
  "main.calling_method": "emptydrops",
  "main.forceCellNum": 0,
  "main.Intron": true
}
```

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
Parameter	Type	Description
main.Outdir	Directory	MANDATORY. Output directory.
main.SampleName	String	MANDATORY. Sample name.
main.cDNA_Fastq1	Fastq	MANDATORY. cDNA Read 1 in fastq format. Can be gzipped. Fastqs from different lanes can be seperated with commas. For example, "L01_read_1.fq.gz,L02_read_1.fq.gz,...".
main.cDNA_Fastq2	Fastq	MANDATORY. cDNA Read 2 in fastq format. Can be gzipped. Fastqs from different lanes can be seperated with commas. For example, "L01_read_2.fq.gz,L02_read_2.fq.gz,...".
main.Oligo_Fastq1	Fastq	MANDATORY. Oligo Read 1 in fastq format. Can be gzipped. Fastqs from different lanes can be seperated with commas. For example, "L01_oligo_1.fq.gz,L02_oligo_1.fq.gz,...".
main.Oligo_Fastq2	Fastq	MANDATORY. Oligo Read 2 in fastq format. Can be gzipped. Fastqs from different lanes can be seperated with commas. For example, "L01_oligo_2.fq.gz,L02_oligo_2.fq.gz,...".
main.BeadsBarcode	json	MANDATORY. cDNA Read structure configure and whitelist file.
main.OligoBarcode	json	MANDATORY. oligo Read structure configure and whitelist file.
main.Root	Directory	MANDATORY. Directory of this pipeline.
main.Refdir	Directory	MANDATORY. STAR index directory of genome reference.
main.Gtf	File Path	MANDATORY. gtf file of genome reference.
main.Oligo_type8	File Path	MANDATORY. Whitelist of oligo.

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改为  
oligomix

```
"main.OligoBarcode": "/Local/path/to/pipeline/DNBC4tools/config/DNBelabC4_scRNA_oligo_readStructure.json",
```

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main.Species	String	Optional, default: NA. Species.
main.expectCellNum	Integer	Optional, default: 3000, expected number of recovered beads for emptydrops.
main.calling_method	String	Optional, default: emptydrops, cell calling method, choose from barcoderanks and emptydrops.
main.forceCellNum	Integer	Optional, default: 0, force pipeline to use this number of beads. 0 means do not use "forceCellNum" to cut off.
main.Intron	Boolean	Optional, default: true, true or false include intronic reads in count.
main.mtgenes	String	Optional, default: auto, set mitochondrial genes (mtgene list file path) or auto. mtgenes's structure like <a href="#">this</a>
main.clusterdim	Integer	Optional, default: 20, the principal components used for clustering.
main.doublepercentage	Float	Optional, default: 0.05, assuming doublet formation rate, tailor for your dataset.
main.mitpercentage	Integer	Optional, default: 15, filter cells with mtgenes percentage.
main.minfeatures	Integer	Optional, default: 200, filter cells with minimum nfeatures.
main.PCusage	Integer	Optional, default: 50, the total number of principal components for PCA.
main.resolution	Float	Optional, default: 0.5, cluster resolution.

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Required parameter:

--name sample name.

--cDNAfastq1 The R1 sequence of the sample cDNA library, multiple samples are separated by commas, and the sequence is the same as that of cDNA R2.

--cDNAfastq2 The R2 sequence of the sample cDNA library, multiple samples are separated by commas, and the sequence is the same as that of cDNA R1.

--oligofastq1 The R1 sequence of the sample oligo library, multiple samples are separated by commas, and the sequence is the same as that of oligo R2.

--oligofastq2 The R2 sequence of the sample oligo library, multiple samples are separated by commas, and the sequence is the same as that of oligo R1.

--starIndexDir The STAR index path of the species' genome.

--gtf: Species annotation file, which needs to match the genome file.

Optional parameter:

--species The species name corresponding to the sample, the default is NA, it is recommended to select, the species name will be displayed in the report, and if the species is human and mouse, the cell population will be annotated in the analysis.

--outdir The result is a directory on which a directory of sample names is generated. Defaults to the current path.

--thread The number of processes used for software analysis, the default is 4.

--cDNAconfig The structure and whitelist of cell barcode file of the cDNA library, default is the config location of the package.

--oligoconfig The structure and whitelist of cell barcode file of the oligo library, default is the config location of the package.

--oligotype The whitelist file of the oligo library, which defaults to the config location of the package.

--calling\_method Methods for identifying empty droplets, including barcoderanks and emptydrops, default is emptydrops.

--expectcells The expected number of beads to capture, this parameter applies to emptydrops, defaults is 3000.

--forcecells Cut off the specified number of beads for subsequent analysis, default is 0, no cut off.

--mtgenes Set mitochondrial genes(mtgene list file path) or auto, auto is find genes starting with MT or mt. mtgenes's file like [this](#)

--process Select the steps to be analyzed, including data, count, analysis, and report. If this parameter is selected, some of the steps are in error and the step can be rerun.

FLAG parameter:

--no\_introns By default, the reads of exon and intron are calculated at the same time. If this parameter is selected, the intron reads will not be calculated. It is recommended not to select.

--mixseq cDNA and oligo library are sequenced on the same chip, and this parameter is added when the sequencing mode is cDNA sequencing mode. In the presence of this parameter, the --oligoconfig parameter is invalid, and

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

### summary&质控图

#### Summary ?

Sample name	PBMC
Species	human
Estimated number of cells	10,988
Mean reads per cell	5,068
Mean UMI count per cell	3,042
Median UMI counts per cell	2,458
Total genes detected	33,132
Mean genes per cell	1,592
Median genes per cell	1,421
Fraction Reads in cells	46.36%
Sequencing saturation	20.47%

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Low quality+failed barcode=total reads (100%) -readspassQC

#### Sequencing ?

##### mRNA

Number of reads	240,126,272
Reads pass QC	61.58%
Reads with exactly matched barcodes	54.21%
Reads with failed barcodes	38.01%
Reads filtered on low quality	0.41%
Q30 bases in cell barcode	92.81%
Q30 bases in UMI	90.46%
Q30 bases in reads	90.50%

##### Droplet index

Number of Reads	105,030,454
Reads pass QC	87.29%
Reads with exactly matched barcodes	77.85%
Reads with failed barcodes	11.26%
Reads filtered on low quality	1.45%
Q30 bases in cell barcode	93.40%
Q30 bases in reads	93.90%

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#### Mapping & Annotation

Raw reads	417,346,234
Mapped reads	393,169,754 (94.21%)
Plus strand	204,275,142
Minus strand	188,894,612
Mitochondria ratio	0.00%
Mapping quality corrected reads	11,111,522

>80%


Reads mapped to genome (Map Quality >= 0)	94.2%
Reads mapped to exonic regions	74.3%
Reads mapped to intronic regions	8.3%
Reads mapped to both exonic and intronic regions	1.1%
Reads mapped antisense to gene	10.4%
Reads mapped to intergenic regions	5.4%
Reads mapped to gene but failed to interpret type	0.5%

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Building No.11, Beishan Industrial Zone, Yantian District, Shenzhen 518083, China

MGI-service@mgi-tech.com | www.mgi-tech.com